

ACCEPTANCE TEST PROCEDURE

For

ANAMORPHIC ATTACHMENT

# 2 (BRASS)

For

High Power Stereoviewer

STAT

Test performed by \_\_\_\_\_ DATE \_\_\_\_\_

\_\_\_\_\_  
DATE \_\_\_\_\_

Anamorphic Attachment Accepted \_\_\_\_\_ DATE \_\_\_\_\_

STAT

Revised 11/7/67

NGA Review Complete

## Acceptance Test Procedure for Anamorphic Attachment


The tests will be performed using the standard  High Power Stereoviewer equipped with the  3X, 6X and 10X objectives, the  1.3X objective and the  6X and 10X eyepieces except as noted.

STAT

STAT

## 1. Resolution

Resolution will be measured axially and at four places at the edge of the field as illustrated in the sketch of the field.



The resolution values of the HPSV without the Anamorphic Attachment will be considered as the reference values. The resolution read with the Anamorphic Attachment in place will be compared with the reference values. The resolution values of the HPSV with Anamorphic Attachment should be at least 90% of the resolution values of the HPSV. A high contrast, black bars on clear background, target will be used.

			HPSV Resolution at Field Position					HPSV With Anamorphic Attachment					Accept	Reject
10X Eyepiece			A	L	R	T	B	A	L	R	T	B		
TARGET BLOCKS 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 125, 150, 175, 200.	1.3X	1X	70	20	20	30	30	80	20	20	30	30		
		2X	125	30	30	40	50	125	30	20	40	30		
	3X	1X	165	54	54	27	27	165	54	54	54	82		
		2X	275	110	137	110	137	275	110	192	110	137		
	6X	1X	343	165	192	137	247	275	165	192	165	275		

BLOCKS IN RESOLUTION TARGET

A L R  $\frac{-2-}{T}$  B

6X	2X	420	180	300	180	300	420	240	300	240	300		
(NOMINAL) REFOCUS EDGE		180	540	480	300	420	180	540	540	600	540		
10X FOCUS CENTER 1X		600	180	180	180	240	600	180	240	240	300		
REFOCUS EDGE EACH TARGET		-	480	300	360	360	-	480	420	420	420		
10X FOCUS CENTER 2X		600	-	-	-	-	600	240	180	240	300		
FOCUS CENTER													
REFOCUS AT EDGE	2X	900	540	600	600	480	900	600	600	600	540		

Comments:

FOR 6X & 10X  
OBS. TARGET BLOCKS EDGE OF CIRCLE OF TARGET AT EDGE OF FIELD.  
60, 120, 180, 240, 300, 360, 420, 480, 540, 600, 750, 900, 1050, 1200

			HPSV Resolution at Field Position					HPSV With Anamorphic Attachement						
6X Eyepiece Objectives	Zoom	Setting	A	L	R	T	B	A	L	R	T	B	Accept	Reject
1.3X		1X	60	20	20	10	20	60	20	20	10	20		
1.3X		2X	100	40	30	20	30	100	40	30	40	50		
3X		1X	137	54	54	54	54	165	110	54	54	82		
3X		2X	247	165	137	110	137	247	165	137	110	165		
6X		1X	275	110	165	110	192	275	165	137	110	192		
6X		2X	412	343	343	343	343	412	343	343	275	343		
10X		1X	420	240	240	240	300	420	300	240	240	300		
10X		2X	750	540	420	360	420	600	540	420	420	480		

Comments:

Focus center & focus

-3-

## 2. Field Size

A scale will be placed in the object plane and the field size will be measured. The Anamorphic Attachment shall not cause more than a 5% loss of field when compared with the standard HPSV.

		HPSV	HPSV With Anamorphic Attachment		
10X Eyepiece Objectives	Zoom Setting	Field Size in mm		Accept	Reject
1.3X	1X	<u>13.3</u>	<u>13.3</u>	_____	_____
1.3X	2X	<u>6.75</u>	<u>6.7</u>	_____	_____
3X	1X	<u>6.5</u>	<u>6.4</u>	_____	_____
3X	2X	<u>3.0</u>	<u>3.05</u>	_____	_____
6X	1X	<u>3.25</u>	<u>3.2</u>	_____	_____
6X	2X	<u>1.52</u>	<u>1.55</u>	_____	_____
10X	1X	<u>1.9</u>	<u>1.87</u>	_____	_____
10X	2X	<u>.9</u>	<u>.93</u>	_____	_____

Comments:

		HPSV	HPSV With Anamorphic Attachment		
Objectives	6X Eyepiece Zoom Setting	Field Size in mm		Accept	Reject
1.3X	1X	<u>12.9</u>	<u>12.7</u>	_____	_____
1.3X	2X	<u>6.55</u>	<u>6.5</u>	_____	_____
3X	1X	<u>6.25</u>	<u>6.1</u>	_____	_____
3X	2X	<u>3.04</u>	<u>2.94</u>	_____	_____
6X	1X	<u>3.15</u>	<u>3.1</u>	_____	_____
6X	2X	<u>1.52</u>	<u>1.5</u>	_____	_____
10X	1X	<u>1.85</u>	<u>1.83</u>	_____	_____

10X Approved For Release 2004/11/30 : CIA-RDP78B04770A000700020020-2  
Comments:

-4-

## 3. Anamorphic Magnification

STAT  
STAT

In this test a ☐ 10X wide field eyepiece will be used instead of the ☐ eyepiece. Its purpose is to accept a scale which will be used for measuring the lengths of perpendicular meridians. A suitable scale or grid will be used in the object plane. The ratio of the lengths of perpendicular meridians is a measure of the Anamorphic Magnification. The Anamorphic Magnification range shall be from 1.0 to 2.2X.

## HPSV WITH ANAMORPHIC ATTACHMENT

		Calculated Anamorphic Magnification (Ratio of Perpendicular Meridians)	Accept	Reject
Anamorphic Scale Setting				
3X obj.	1.0			
	1.2	$\frac{10}{8} = 1.25$		
1X Zoom Setting	1.4	$\frac{10}{7} = 1.43$		
	1.6	$\frac{10}{6.3} = 1.588$		
	1.8	$\frac{10}{5.5} = 1.82$		
	2.0	$\frac{10}{4.9} = 2.04$		
	2.2	$\frac{10}{4.5-4.6} = 2.17 - 2.22$		

Comments:

-5-

## 4. Eye Point Extension and Eye Relief

10X  EYEPiece

STAT

The difference in length between the standard HPSV eyepoint and the eyepoint of the HPSV with Anamorphic Attachment will be calculated.

Both measurements will be made relative to a fixed point on the HPSV.

		Accept	Reject
Distance with Anamorphic Attachment	<u>177 mm</u>	<u>          </u>	<u>          </u>
Distance with Standard HPSV	<u>153 mm</u>	<u>          </u>	<u>          </u>
Difference - Eyepoint Extension	<u>24 mm</u>	<u>          </u>	<u>          </u>

The eyepoint extension shall be no more than 3 inches.

The eye relief shall be measured from the exit pupil to the eyepiece.

STAT

10X  EYEPiece

	Standard HPSV	HPSV with Anamorphic Attachment	Accept	Reject
Eye Relief	<u>16</u>	<u>15.5 mm</u>	<u>          </u>	<u>          </u>
6X <input type="text"/> EYEPiece	<u>21 mm</u>			

STAT

The Interpupillary separation shall be measured with and without the Anamorphic Attachments in place.

	Standard HPSV	HPSV with Anamorphic Attachment	Accept	Reject
IPD	<u>54.5 - 73</u>	<u>53.5 - 71.5</u>	<u>          </u>	<u>          </u>

Comments:

-6-

## 5. Interchangeability

The time required to remove the Anamorphic Attachment shall be less than five minutes, without the use of special tools.

Time Required for Removal of the Anamorphic Attachment

Accept      Reject

30 SECONDS ~~Minutes~~

Comments:

## 6. Anamorphic Axis Orientation

Verification will be made that the direction of anamorphic magnification shall be rotatable through 360°.

Accept      Reject

X

Comments:

## 7. Percent Transmission

The transmission of the Anamorphic Attachment shall be determined.

A small diameter collimated beam of light will be transmitted through the Anamorphic Attachment equipped with the ☐ 10X wide field compensating eyepiece. The light energy will be measured and will be compared to the light energy passing through the ☐ 10X eyepiece. The ratio of the two values obtained will be a measure of the light

STAT

STAT

-7-

transmission of the Anamorphic Attachment.

	Light	Energy
	1X	2.2X
(1) Anamorphic Attachment with 10X	<div style="border: 1px solid black; width: 40px; height: 40px; display: inline-block;"></div>	Eyepiece
(2) 10X		

STAT

$$\% \text{ Transmission} = \frac{(1)}{(2)} \times 100$$

=

90%8. PERMANENT RUNOUT MEASURED IN EYEPIECE FOCAL PLANE

HPSU WITH ANAMORPHIC ATTACHMENT

HPSU

1.2mm

0.2mm

9. NO PERCEPTIBLE CHANGE WHEN ROTATING PLANE OF ANAMORPHICA AT 1X SETTING

10. REMAINS IN FOCUS THROUGHOUT MICROSCOPE ZOOM RANGE &amp; ANAMORPHIC ZOOM RANGE



ACCEPTANCE TEST PROCEDURE

For

ANAMORPHIC ATTACHMENT #1

For

High Power Stereoviewer



STAT

Test performed by \_\_\_\_\_ DATE \_\_\_\_\_

\_\_\_\_\_  
DATE \_\_\_\_\_

Anamorphic Attachment Accepted \_\_\_\_\_ DATE \_\_\_\_\_



STAT


Revised 11/7/67

# Acceptance Test Procedure for Anamorphic Attachment

The tests will be performed using the standard  High Power Stereoviewer equipped with the  5X, 6X and 10X objectives, the  1.3X objective and the  6X and 10X eyepieces except as noted.

## 1. Resolution

Resolution will be measured axially and at four places at the edge of the field as illustrated in the sketch of the field.



The resolution values of the HPSV without the Anamorphic Attachment will be considered as the reference values. The resolution read with the Anamorphic Attachment in place will be compared with the reference values. The resolution values of the HPSV with Anamorphic Attachment should be at least 90% of the resolution values of the HPSV. A high contrast, black bars on clear background, target will be used.

		HPSV Resolution at Field Position					HPSV With Anamorphic Attachment					Accept	Reject
Objectives	10X Eyepiece Zoom Setting	A	L	R	T	B	A	L	R	T	B		
1.3X	1X	-	-	-	-	-	90	20	10	10	20		
1.3X	2X	-	-	-	-	-	125	40	40	40	30		
3X	1X	-	-	-	-	-	192	110	54	54	82		
3X	2X	-	-	-	-	-	275	165	137	110	137		
6X	1X	-	-	-	-	-	360	180	180	180	180		

See data on

-2-

6X	2X	—	—	—	—	—	480	240	300	240	300	—	—
10X	1X	—	—	—	—	—	600	480	240	300	420	—	—
10X	2X	—	—	—	—	—	750	600	600	600	600	—	—

See data on #2

Comments:

		HPSV Resolution at Field Position					HPSV With Anamorphic Attachement					Accept	Reject
6X Eyepiece Objectives	Zoom Setting	A	L	R	T	B	A	L	R	T	B		
1.3X	1X	—	—	—	—	—	50	20	10	10	20	—	—
1.3X	2X	—	—	—	—	—	90	40	40	40	40	—	—
3X	1X	—	—	—	—	—	137	82	54	54	82	—	—
3X	2X	—	—	—	—	—	247	165	137	110	137	—	—
6X	1X	—	—	—	—	—	225	165	165	165	192	—	—
6X	2X	—	—	—	—	—	412	275	343	247	343	—	—
10X	1X	—	—	—	—	—	480	360	300	240	300	—	—
10X	2X	—	—	—	—	—	750	600	540	600	600	—	—

See data on #2

Comments:

-3-

## 2. Field Size

A scale will be placed in the object plane and the field size will be measured. The Anamorphic Attachment shall not cause more than a 5% loss of field when compared with the standard HPSV.

		HPSV	HPSV With Anamorphic Attachment		
10X Eyepiece		Field Size in mm		Accept	Reject
Objectives	Zoom Setting				
1.3X	1X	_____	13.3 13.5	_____	_____
1.3X	2X	_____	6.7 6.36	_____	_____
3X	1X	_____	6.35 6.1	_____	_____
3X	2X	_____	3.07 2.95	_____	_____
6X	1X	_____	3.2 3.05	_____	_____
6X	2X	_____	1.55 1.48	_____	_____
10X	1X	_____	1.87 1.78	_____	_____
10X	2X	_____	.9 .89	_____	_____

Comments:

*See data on #2**10X Wild  
E.P.**a different  
10X Wild E.P.*

		HPSV	HPSV With Anamorphic Attachment		
6X Eyepiece		Field Size in mm		Accept	Reject
Objectives	Zoom Setting				
1.3X	1X	_____	13.2 13.0	_____	_____
1.3X	2X	_____	6.32 6.5	_____	_____
3X	1X	_____	6.0 6.12	_____	_____
3X	2X	_____	2.9 2.95	_____	_____
6X	1X	_____	3.03 3.1	_____	_____
6X	2X	_____	1.47 1.5	_____	_____
10X	1X	_____	1.79 1.82	_____	_____

10X Approved For Release 2004/11/30 : CIA-RDP78B04770A000700020020-2  
Comments:

*See data on #2**Second  
6X E.P.**6X E.P. used before on #2*

-4-

## 3. Anamorphic Magnification

STAT In this test a  10X wide field eyepiece will be used instead of the  
STAT  eyepiece. Its purpose is to accept a scale which will be used  
for measuring the lengths of perpendicular meridians. A suitable scale  
or grid will be used in the object plane. The ratio of the lengths of  
perpendicular meridians is a measure of the Anamorphic Magnification.  
The Anamorphic Magnification range shall be from 1.0 to 2.2X.

## HPSV WITH ANAMORPHIC ATTACHMENT

		Calculated Anamorphic Magnification (Ratio of Perpendicular Meridians)	Accept	Reject
Anamorphic Scale Setting				
3X obj.	1.0			
	1.2	$\frac{10}{8.2} = 1.22$		
1X Zoom Setting	1.4	$\frac{10}{7.15} = 1.40$		
	1.6	$\frac{10}{6.3} = 1.587$		
	1.8	$\frac{10}{5.6} = 1.786$		
	2.0	$\frac{10}{5.0} = 2.00 +$		
	2.2	$\frac{10}{4.5} = 2.22$		

Comments:

-5-

## 4. Eye Point Extension and Eye Relief

The difference in length between the standard HPSV eyepoint and the eyepoint of the HPSV with Anamorphic Attachment will be calculated. Both measurements will be made relative to a fixed point on the HPSV.

		Accept	Reject
Distance with Anamorphic Attachment	_____	_____	_____
Distance with Standard HPSV	_____	_____	_____
Difference - Eyepoint Extension	_____	_____	_____

The eyepoint extension shall be no more than 3 inches.

The eye relief shall be measured from the exit pupil to the eyepiece.

	Standard HPSV	HPSV with Anamorphic Attachment	Accept	Reject
Eye Relief	_____	_____	_____	_____

The Interpupillary separation shall be measured with and without the Anamorphic Attachments in place.

	Standard HPSV	HPSV with Anamorphic Attachment	Accept	Reject
IPD	_____	_____	_____	_____

Comments:

-6-

# 5. Interchangeability

The time required to remove the Anamorphic Attachment shall be less than five minutes, without the use of special tools.

Time Required for Removal of the Anamorphic Attachment

Accept Reject

Minutes

Comments:

# 6. Anamorphic Axis Orientation

Verification will be made that the direction of anamorphic magnification shall be rotatable through 360°.

Accept Reject

X

Comments:

# 7. Percent Transmission

The transmission of the Anamorphic Attachment shall be determined.

A small diameter collimated beam of light will be transmitted through the Anamorphic Attachment equipped with the ☐ 10X wide field compensating eyepiece. The light energy will be measured and will be compared to the light energy passing through the ☐ 10X eyepiece. The ratio of the two values obtained will be a measure of the light

STAT

STAT

-7-

transmission of the Anamorphic Attachment.

Light	Energy
1X	2.2X

- |     |                                |  |   |
|-----|--------------------------------|--|---|
| (1) | Anamorphic Attachment with 10X | <div style="border: 1px solid black; width: 40px; height: 40px; display: inline-block;"></div> | Eyepiece  |
| (2) |                                | 10X  | <div style="border: 1px solid black; width: 40px; height: 40px; display: inline-block;"></div> Eyepiece |

STAT

$$\% \text{ Transmission} = \frac{(1)}{(2)} \times 100 = \underline{\underline{90\%}}$$

9. PECHAN RUN OUT MEASURED IN EYEPIECE FOCAL PLANE  
 HPS WITH ANAMORPHIC ATTACHMENT HPS  
 0.4 mm 0.2 mm

9. NO PERCEPTIBLE CHANGE WHEN ROTATING PLANE OF  
 ANAMORPHISM AT 1X SETTING.

10. REMAINS IN FOCUS THROUGHOUT MICROSCOPE FOCUS RANGE  
 & ANAMORPHIC ZOOM RANGE